

BRODA[®]

Vanguard

Positioning Wheelchair



Shown with the Acute Care Configuration (ACC) full vinyl wrap option

Patient, Facility, & Caregiver Benefits

- Pressure redistribution that aids in maintaining skin integrity
- Improved air flow, decreasing heat and moisture build up
- Accommodation for conditions such as kyphosis and lordosis



Anterior and posterior tilt aids in safe patient handling



Swing Away, Removable Arm Supports facilitate side-transfers



Manual foot pump allows safe, no lift tilt positioning




Concave back for gluteal shelf improves user support & dignity

brodaseating.com | 1-844-552-7632

BARIATRIC POSITIONING
CARE REDEFINED

UPHOLSTERY OPTIONS



FRAME OPTIONS



SPECIFICATIONS

Degrees of Tilt	up to 10° Anterior, up to 30° Posterior
Degrees of Recline	up to 45°
Weight Capacity	up to 600 lbs. (272 kg)
Seat Widths	24", 26", 28", 30" (61cm, 66cm, 71cm, 76cm)
Seat Depths	21", 23" (53cm, 58cm)
Seat Heights	20", 21", 22", 23", 24" (50cm, 53cm, 55cm, 58cm, 60cm)
Wheels	6" Front Casters: Total Lock (15cm) 8" Rear Casters: Fixed Multi-Directional Locking (20cm)
Lifetime Frame Warranty	10 years on 16 Gauge tubular steel frame
Pressure-Reducing Seat & Backrest	Integrated with Comfort Tension Seating®
Standard Padding	Headrest cushion, upper and lower side pads, shoulder bolsters, seat and back pad with cover, calf pad

OPTIONAL ACCESSORIES



O² Holder



Padded Thigh Strap



Half Tray