

Midline

Positioning Wheelchair



Shown with rear-mounted mag wheel option

Patient, Facility, & Caregiver Benefits

- Increases quality of life and engagement with surroundings
- Allows for easy maneuverability and safe patient adjustments
- Reduces risk of injury and falls with supportive stable positioning



Full tilt and recline functionality provides superior positioning



Adjustable seat surface height adapts the Midline to support users at multiple care levels



Tilt-in-space enhances comfort



Comfort Tension Seating® increases support as well as sitting tolerance



Midline

Positioning Wheelchair

ENHANCED POSITIONING WITH REFINED FEATURES

UPHOLSTERY OPTIONS







FRAME OPTIONS





OPTIONAL ACCESSORIES



Foot Box



O² Holder



Lateral Support



Half Tray

SPECIFICATIONS

| Degrees of Tilt | Up to 38° |
|--------------------------------------|--|
| Recline | Up to 90° (Trendelenburg capable) |
| Weight Capacity | up to 350 lbs. (114 kg) |
| Seat Widths | 16", 18", 20", 22", 24", 26", 28" (40cm, 45cm, 50cm, 55cm, 61cm, 66cm, 71cm) |
| Seat Depths | 15.5", 17", 18.5", 20" (39cm, 43cm, 45cm, 58cm) |
| Floor-to-Seat Heights | Adjustable 17", 18", 19", 20", 21", 22" (43cm - 55cm) |
| Wheels | 5" Front Casters: Total Lock (12.5cm) 5" Rear Casters: Directional or Total Lock (20cm) |
| Lifetime Warranty Frame | 10 years on 16 Gauge tubular steel frame |
| Pressure-Reducing Seat & Backrest | Integrated with Comfort Tension Seating® |
| Standard Padding | Neckrest, shoulder bolsters, upper and lower side pads, seat and back pad, calf and sole pad |







Padded Thigh Strap